

EQUIVET Ocular Lavage Kit

The subpalpebral ocular lavage catheter is designed to keep a lavage catheter in situ for an extended period of time with minimum discomfort to the patient and operator. Medication can readily be administered via the long tubing. The proximal tip of the catheter is soft and designed to fit comfortably against the conjunctiva. With the needles injection site it is a lot easier for the client to treat the horse. If using a thicker medication choose the set with the 8 Fr size silicon catheter instead of the 5 Fr size silicon catheter. The design and length of the trocar makes it easier and safer to place the catheter.



EQUIVET



EQUIVET Ocular Lavage 5 Fr Kit

Components:


Silicon Lavage Catheter 5 Fr with 0,76 mm I.D. and 1,7 mm O.D., 90 cm long

Trocar for 5 Fr

Needle 21G x 1" blunt

Needleless Injection Site with protection cap

Cat. no. 230700



EQUIVET Ocular Lavage 8 Fr Kit

Components:

Silicon Lavage Catheter 8 Fr with 1,5 mm I.D. and 2,5 mm O.D., 90 cm long

Trocar for 8 Fr

Needle 16G x 1" blunt


Needleless Injection Site with protection cap

Cat. no. 230705



For further instruction about

"How to use the EQUIVET Ocular Sub-Palpebral Lavage set"
see our website at www.kruuse.com.



EQUIVET

EXPERIENCE · QUALITY · SERVICE
www.kruuse.com

